

Carbon Pin Type Mount

	<p>型番: 01608X-BA(10) 01608S-BA(10)</p>	<p>Table diameter: 14.3 mm Pin diameter: 6.4 mm Pin length: 16 mm Surface: Standard, Ultra-Smooth finish Compatible SEM: AMRAY 900</p>
	<p>型番: 01606X-BA(10) 01606X-MB(100)</p>	<p>Table diameter: 12.7 mm Pin diameter: 3.2 mm Pin length: 8 mm Surface: Standard finish Compatible SEM: LEO/Cambridge/Philips/ETEC/FEI</p>
	<p>型番: 01609X-BA(10) 01609X-BA(10)</p>	<p>Table diameter: 12.7 mm Pin diameter: 3.2 mm Pin length: 15 mm Surface: Standard finish Compatible SEM: AMRAY</p>