



















Round Cylindrical Mount

	<p>型番: 1538L-BA(10), 1538L-MB(100) サイズ: 32 mm Dia x 10 mm High</p>	<p>Diameter: 32 mm Height: 10 mm Surface:</p> <p>Compatible SEM: JEOL</p>
	<p>型番: 01501-BA(10), 01501-MB(100) 型番: 01501L-BA(10), 01501L-MB(100) 型番: 01501P-BA(10), 01501P-MB(100) サイズ: 9.5 mm Dia x 9.5 mm High</p>	<p>Diameter: 9.5 mm Height: 9.5 mm Surface: Lathe finish</p> <p>Compatible SEM: JEOL</p>
	<p>型番: 01502-BA(10), 01502-MB(100) 型番: 01502L-BA(10), 01502L-MB(100) 型番: 01502P-BA(10), 01502P-MB(100) サイズ: 12 mm Dia x 5 mm, High</p>	<p>Diameter: 12 mm Height: 5 mm Surface: Lathe finish</p> <p>Compatible SEM: JEOL 840 SERIES (Including Models 5800 & 6400)</p>
	<p>型番: 01503-BA(10), 01503-MB(100) 型番: 01503L-BA(10), 01503L-MB(100) 型番: 01503P-BA(10), 01503P-MB(100) サイズ: 12 mm Dia x 10 mm High</p>	<p>Diameter: 12 mm Height: 10 mm Surface: Lathe finish</p> <p>Compatible SEM: JEOL 840 SERIES (Including Models 5800 & 6400)</p>

 <p>Technical drawing of a cylindrical magnet with a diameter of 15.00 mm and a height of 10.00 mm. The drawing includes a top view, a side view, and a perspective view.</p>	 <p>型番: 01504-BA(10), 01504-MB(100) 型番: 01504L-BA(10), 01504L-MB(100) 型番: 01504P-BA(10), 01504P-MB(100) サイズ: 15 mm Dia x 10 mm, High</p>	<p>Diameter: 15 mm Height: 10 mm Surface: Lathe finish</p> <p>Compatible SEM: TOPCON/ISI</p>
 <p>Technical drawing of a cylindrical magnet with a diameter of 25.40 mm and a height of 20.00 mm. The drawing includes a top view, a side view, and a perspective view.</p>	 <p>型番: 01510-BA(10), 01510-MB(100) 型番: 01510L-BA(10), 01510L-MB(100) サイズ: 25 mm Dia x 20 mm, High</p>	<p>Diameter: 25 mm Height: 20 mm Surface: Lathe finish</p> <p>Compatible SEM: 1" Metallographic Universal</p>
 <p>Technical drawing of a cylindrical magnet with a diameter of 25.4 mm and a height of 13 mm. The drawing includes a top view, a side view, and a perspective view.</p>	 <p>型番: 01511-BA(10), 01511-MB(100) 型番: 01511L-BA(10), 01511L-MB(100) サイズ: 25 mm Dia x 13 mm, High</p>	<p>Diameter: 25 mm Height: 13 mm Surface: Lathe finish</p> <p>Compatible SEM: 1" Metallographic Universal</p>
 <p>Technical drawing of a cylindrical magnet with a diameter of 25.4 mm and a height of 8 mm. The drawing includes a top view, a side view, and a perspective view.</p>	 <p>型番: 01512-BA(10), 01512-MB(100) 型番: 01512L-BA(10), 01512L-MB(100) サイズ: 25 mm Dia x 8 mm, High</p>	<p>Diameter: 25 mm Height: 8 mm Surface: Lathe finish</p> <p>Compatible SEM: 1" Metallographic Universal</p>
 <p>Technical drawing of a cylindrical magnet with a diameter of 15.00 mm and a height of 5.00 mm. The drawing includes a top view, a side view, and a perspective view.</p>	 <p>型番: 01514-BA(10), 01514-MB(100) 型番: 01514L-BA(10), 01514L-MB(100) 型番: 01514P-BA(10), 01514P-MB(100) サイズ: 15 mm Dia x 5 mm, High</p>	<p>Diameter: 15 mm Height: 5 mm Surface: Lathe finish</p> <p>Compatible SEM: TOPCON/ISI</p>

<p> ϕ 15.00 mm 15.00 mm </p>	<p> 型番: 01515-BA(10), 01515-MB(100) 型番: 01515L-BA(10), 01515L-MB(100) 型番: 01515P-BA(10), 01515P-MB(100) サイズ: 15 mm Dia x 15 mm, High </p>	Diameter: 15 mm Height: 15 mm Surface: Lathe finish Compatible SEM: TOPCON/ISI
<p> ϕ 10 mm 5 mm </p>	<p> 型番: 01516-BA(10), 01516-MB(100) サイズ: 10 mm Dia x 5 mm, High </p>	Diameter: 10 mm Height: 5 mm Surface: Lathe finish Compatible SEM: JEOL
<p> ϕ 32.00 mm 22.50 mm </p>	<p> 型番: 01536L-BA(10) サイズ: 32 mm Dia x 22 mm, High </p>	Diameter: 32 mm Height: 22 mm Surface: Lathe finish Compatible SEM: JEOL
<p> ϕ 10 mm 5 mm </p>	<p> 型番: 01537-BA(10), 01537-MB(100) 型番: 01537L-BA(10), 01537L-MB(100) サイズ: 10 mm Dia x 5 mm, High </p>	Diameter: 10 mm Height: 5 mm Surface: Lathe finish Compatible SEM: JEOL 840 SERIES
<p> ϕ 32 mm 10 mm </p>	<p> 型番: 01538-BA(10), 01538-MB(100) 型番: 01538L-BA(10), 01538L-MB(100) サイズ: 32 mm Dia x 10 mm, High </p>	Diameter: 32 mm Height: 10 mm Surface: Lathe finish Compatible SEM: JEOL

	 <p>型番: 01539-BA(10), 01539-MB(100) 型番: 01539L-BA(10), 01539L-MB(100) サイズ: 25 mm Dia x 5 mm, High</p>	<p>Diameter: 25 mm Height: 5 mm Surface: Lathe finish</p> <p>Compatible SEM: JEOL</p>
	 <p>型番: 01577-BA(10), 01577-MB(100) 型番: 01577L-BA(10), 01577L-MB(100) サイズ: 26 mm Dia x 10 mm, High</p>	<p>Diameter: 26 mm Height: 10 mm Surface: Lathe finish</p> <p>Compatible SEM: HITACHI</p>
	 <p>型番: 01578-BA(10), 01578-MB(100) 型番: 01578L-BA(10), 01578L-MB(100) サイズ: 14 mm Dia x 10 mm, High</p>	<p>Diameter: 14 mm Height: 10 mm Surface: Lathe finish</p> <p>Compatible SEM: HITACHI</p>
	 <p>型番: 01579-BA(10), 01579-MB(100) 型番: 01579L-BA(10), 01579L-MB(100) サイズ: 15 mm Dia x 13 mm, High</p>	<p>Diameter: 15 mm Height: 13 mm Surface: Lathe finish</p> <p>Compatible SEM: HITACHI</p>